

EVIS LUCERA E L I T E

EVIS LUCERA ELITE GASTROINTESTINAL VIDEOSCOPE

GIF-HQ290

Amazing image quality for upper GI endoscopy



OLYMPUS GIF-HQ290










Main Features

Dual Focus

Dual focus, a unique OLYMPUS optical innovation, allows the user to select from two focus settings. With the simple push of a scope button, the desired depth of field for observation can be optimized to either the near field or normal field.






....

Enhanced Image Quality

A new level of resolving power is achieved thanks to a new optical system. Combination of the scope and the system provides superior imaging.

NBI (Narrow Band Imaging)

NBI in EVIS LUCERA ELITE provides twice the viewable distance of EVIS LUCERA SPECTRUM. The greatly improved performance of NBI opens up exciting new clinical applications and reinforces NBI's position as the standard of care for GI endoscopy.





Specifications

Optical System	Field of view	Normal focus mode 140°
	Field of view	Near focus mode 140°
	Direction of view	Forward viewing
	Depth of field	Normal focus mode 7-100 mm
		Near focus mode 3-7 mm
Insertion Section	Distal end outer diameter	10.2 mm
	Objective Lens Air/Water Nozzle Right Auxiliary Water Channel Distal end enlarged Up Light-guide Lens Left Instrument Channel Outlet	
	Down	
	Insertion tube outer diameter 9.9 mm	
	Working length	1030 mm

Water Jet

A forward water jet helps to keep the mucosa clear during observation and treatment in the upper GI tract.



Waterproof One-touch Connector

A new connector design minimizes the effort required for setup prior to and in between cases. In addition, it is fully submersible and eliminates the need for a water-resistant cap and the associated risk of an expensive repair due to accidental immersion.





Channel inner diameter	2.8 mm
Minimum visible distance	3.0 mm from the distal end
Direction from which endotherapy accessories enter and exit the endoscopic image	
Angulation range	Up 210°
	Down 90°
	Right 100°
	Left 100°
1350 mm	
Video System Center OLYMPUS CV-290 Xenon Light Source OLYMPUS CLV-290SL, CLV-290	
	Minimum visible distance Direction from which endotherapy accessories enter and exit the endoscop image Angulation range 1350 mm Video System Center OLYM

Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.

